

New ALPRO-product

**PlastiSept eco**

We have introduced our innovation PlastiSept *eco* at the fair IDS in Cologne. The special feature of this new surface disinfection is that we refrain entirely from using QAC (Quaternary Ammonium Compounds), alcohol, phenols and aldehydes in its formulation. Therefore PlastiSept *eco* is particularly gentle and causes no streaks or discolorations on the disinfected surfaces.

The modern combination of sec. and tert. alcyamines is

- Free of substances damaging the material
- Free of allergenic and sensitizing substances
- Particularly gentle on all materials
- Free of alcohol and QAC
- Effective: Bactericidal incl. TBC and MRSA, yeasticidal and virucidal against enveloped viruses in Carrier Test DVV 2012 incl. HBV, HCV, HIV, herpes simplex, influenza (incl. H1N1, H5N1), BVDV, vaccinia, corona viruses (e. g. MERS-CoV) as well as non-enveloped noro viruses



PlastiSept *eco* will be available in April 2015 in foam bottles, canisters and cubitainers as well as in the form of practical wet wipes in the standard or jumbo size. Of course, PlastiSept *eco* can also be used with the MaxiWipes.

**More information**

AlproZyme and BIB forte *eco*

**High level virucidal processing of instruments**

The combination of AlproZyme and BIB forte *eco* for the processing of instruments is high level virucidal according to DVV 2012. The alkaline-enzymatic pre-cleaning AlproZyme and the following disinfection with BIB forte *eco* are



- extremely economical (each 0.5 % solution)
- very gentle on materials
- VAH certified
- High level virucidal according to DVV 2012

**To our website**

AlproZyme

**New efficacy**

According to RKI (Robert-Koch-Institute) it is recommended to carry out a non protein fixing pre-cleaning of instruments. AlproZyme is especially suitable for this as it has a self-acting high cleaning performance and is virucidal according to DVV 2012. Therefore AlproZyme provides high protection for your staff already when pre-cleaning.

**To AlproZyme**

Update

**New medical program**

Our updated medical program is now available for you. Please feel free to download it on our website.

**Download**

